


Product Name: Recombinant Human FLCN Protein
Catalog #: PEH8409



Summary

Name	FLCN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NFG4
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EKLLGAPTE DTLVQMEKLA DLEEESESWD NSEAEKEEKA PVLPESTEGR ELTQGPAESS SLSGCGSWQP RKLVPFKSLR HMRQVLGAPS FRMLAWHVLM GNQVIWKS RD VDLVQSAFEV LRTMLPVGCV RIIPYSSQYE EAYRCNFLGL SPHVQIPPHV LSSEFAVIVE VHAAARSTLH PVGCEDDQSL SKYEFVVTSG SPVAADRVGP TILNKIEAAL TNQNLSVDVV DQCLVCLKEE WMNKVKVLFK FTKVDSPKE DTQKLLSILG ASEEDNVKLL KFWMTGLSKT YKSHLMSTVR

image



Background

Alternative Names Folliculin (BHD skin lesion fibrofolliculoma protein) (Birt-Hogg-Dube syndrome protein)

Product Name: Recombinant Human FLCN Protein
Catalog #: PEH8409



Note

For Research Use Only , Not for Diagnostic Use