


Product Name: Recombinant Human DMXL2 Protein
Catalog #: PEH8410



Summary

Name	DMXL2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TDJ6
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KNFVDKINLI ERKLFETTAN AHFKVGCPVL ALEVLSKIPK VTKTSALSAK KDQPDFISHR MDDVPSHKA LSDGNGSSGI EWSNVTSSQY DWSQPIVKVD EEPLNLDWGE DHDSALDEEE DDAVGLVMKS TDAREKDKQS DQKASDPNML LTPQEEDDPE GDTEVDVIAE QLKFRACLKI LMTELRTLAT GYEVDGGKLR FQLYNWLEKE IAALHEICNH ESVIKEYSSK TYSKVESDLL DQEEMVDKPD IGSYERHQIE RRRLQAKREH

image



Background

Alternative Names DmX-like protein 2 (Rabconnectin-3)

Product Name: Recombinant Human DMXL2 Protein
Catalog #: PEH8410



Note

For Research Use Only , Not for Diagnostic Use