


**Product Name: Recombinant Human MYLIP Protein**  
**Catalog #: PEH8411**



## Summary

|                              |                                                                                                                                                                                                                                                                                                                                     |
|------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | MYLIP_HUMAN                                                                                                                                                                                                                                                                                                                         |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                                                                 |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                                                                                                                                    |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                                                                                                                                               |
| <b>Accession #</b>           | Q8WY64                                                                                                                                                                                                                                                                                                                              |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                                                                                                                                              |
| <b>Species</b>               | Human                                                                                                                                                                                                                                                                                                                               |
| <b>Molecular Mass</b>        | 34 kDa                                                                                                                                                                                                                                                                                                                              |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                                                                                                                                           |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                                                                                                  |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                                                                                                      |
| <b>Sequence</b>              | VLMEVEVEAK ANGEDCLNQV CRRLGII EVD YFGLQFTGSK GESLWLNLRN<br>RISQQMDGLA PYRLKLRVKF FVEPHLILQE QTRHIFFLHI KEALLAGHLL<br>CSPEQAVELS ALLAQTKFGD YNQNTAKYNY EELCAKELSS ATLNSIVAKH<br>KELEGTSQAS AEYQVLQIVS AMENYGIEWH SVRDSEGQKL LIGVGPEGIS<br>ICKDDFSPIN RIAYPVVQMA TQSGKNVYLT VTKESGNSIV LLFKMISTRA<br>ASGLYRAITE THAFYRCDTV TSAVMMQYSR |

## image



## Background

|                          |                                                                                                                                                                                                      |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | E3 ubiquitin-protein ligase MYLIP (EC 2.3.2.27) (Inducible degrader of the LDL-receptor) (Idol) (Myosin regulatory light chain interacting protein) (MIR) (RING-type E3 ubiquitin transferase MYLIP) |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**Product Name: Recombinant Human MYLIP Protein**  
**Catalog #: PEH8411**

---



**Note**

For Research Use Only , Not for Diagnostic Use