


Product Name: Recombinant Human PIK3R5 Protein
Catalog #: PEH8412



Summary

Name	PI3R5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WYR1
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GEKAGFPGVL DTAKPGKLHT IIPVARYCYT YSWSQDSFDI LQEILLKEQE LLQPGILGDD EEEEEEEEV EEDLETGHC AERDSSLSTS SLASHDSTLS LASSQASGPA LSRHLLTSFV SGLSDGMDSG

image



Background

Alternative Names	Phosphoinositide 3-kinase regulatory subunit 5 (PI3-kinase regulatory subunit 5) (PI3-kinase p101 subunit) (Phosphatidylinositol 4,5-bisphosphate 3-kinase regulatory subunit) (PtdIns-3-kinase regulatory subunit) (Protein FOAP-2) (PtdIns-3-kinase p101) (p101-PI3K)
--------------------------	---

Note

Product Name: Recombinant Human PIK3R5 Protein
Catalog #: PEH8412



For Research Use Only , Not for Diagnostic Use