

Product Name: Recombinant Human TAF15 Protein
Catalog #: PEH8413



Summary

Name	RBP56_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92804
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SQRENYSHHT QDDRRDVSRY GEDNRGYGGS QGGGRGRGGY DKDGRGPMTG SSGGDRGGFK NFGGHRDYGP RTDADSESDN SDNNTIFVQG LGEGVSTDQV GEFFKQIGII KTNKKTGKPM INLYTDKDTG KPKGEATVSF DDPPSAKAAI DWFDGKEFHG NIIKVSFATR RPEFMRGGGS GGGRRGRGGY RGRGGFQGRG GDPKSGDWVC PNPSCGNMNF ARRNSCNQCN EPRPDSRPS GGDFRGRGYG GERGYRGRGG RGGDRGGYGG DRSGGGYGGD RSSGGGYSGD RSGGGYGGDR SGGGYGGDRG

image



Background

Alternative Names	TATA-binding protein-associated factor 2N (68 kDa TATA-binding protein-associated factor) (TAF(II)68) (TAFII68) (RNA-binding protein 56)
--------------------------	--

Product Name: Recombinant Human TAF15 Protein
Catalog #: PEH8413



Note

For Research Use Only , Not for Diagnostic Use