


Product Name: Recombinant Human NEDD4L Protein
Catalog #: PEH8423



Summary

Name	NED4L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96PU5
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WEERKDAKGR TYYVNHNNRT TTWTRPIMQL AEDGASGSAT NSNNHLIEPQ IRRPRSLSSP TVTSLAPLEG AKDSPVRRRAV KDTLNPNQSP QPSPYNPKP QHKVTQSFLP PGWEMRIAPN GRPFFIDHNT KTTTWEDPRL KFPVHMRSKT SLNPNDLGPL PPGWEERLHL DGRTFYIDHN SKITQWEDPR LQNPAITGPA VPYSREFKQK

image



Background

Alternative Names E3 ubiquitin-protein ligase NEDD4-like (EC 2.3.2.26) (HECT-type E3 ubiquitin transferase NED4L) (NEDD4.2) (Nedd4-2)

Note

Product Name: Recombinant Human NEDD4L Protein
Catalog #: PEH8423



For Research Use Only , Not for Diagnostic Use