


Product Name: Recombinant Human TGS1 Protein
Catalog #: PEH8425



Summary

Name	TGS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96RS0
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EMFLFIEERE DCKILCLCSR AFVEDRKLYN LGLKGYIIRD SGNNSGDQAT EEEEGGYSCG TAESHDSKGI GLDESELDSE AELMRSMGLP LQFGRITAHK DFEVSMNTRN KVKIKKKKHQ KKYLDIVQE SWRKEYEEDD ILASDDPSSI EQYENTRTYE LQSKKDTETE NPPVENTLSP

image



Background

Alternative Names Trimethylguanosine synthase (EC 2.1.1.-) (CLL-associated antigen KW-2) (Cap-specific guanine-N2 methyltransferase) (Hepatocellular carcinoma-associated antigen 137) (Nuclear receptor coactivator 6-interacting protein) (PRIP-interacting protein with methyltransferase motif) (PIMT) (PIPMT)

Product Name: Recombinant Human TGS1 Protein
Catalog #: PEH8425



Note

For Research Use Only , Not for Diagnostic Use