


Summary

Name	ATX2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99700
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LEPWDAGELT ANEELEALEN DVSNGWDPND MFRYNEENYG VVSTYDSSLS SYTVPLERDN SEEFKREAR ANQLAEEIES SAQYKARVAL ENDDRSEEEK YTAVQRNSSE REGHSINTRE NKYIPPGQRN REVISWGSGR QNSPRMGQPG SGSMPSRSTS HTSDFNPNSG SDQRVVNGGV PWPSPCSPS SRPPSRYSQSG PNSLPPRAAT PTRPPSRPPS RPSRPPSHPS AHGSPAPVST MPKRMSSEGP PRMSPKAQRH PRNHRVSAGR GSISGLEFV SHNPPSEAAAT

image



Background

Alternative Names	Ataxin-2 (Spinocerebellar ataxia type 2 protein) (Trinucleotide repeat-containing gene 13 protein)
--------------------------	--

Product Name: Recombinant Human ATXN2 Protein
Catalog #: PEH8427



Note

For Research Use Only , Not for Diagnostic Use