


Summary

Name	FANCJ_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BX63
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NIQDNESTLE VTSLKYSTSP YLLEAASHLS PENFVEDEAK ICVQELQCPK IITKNSPLPS SIISRKEKND PVFLEEAGKA EKIVISRSTS PTFNKQTKRV SWSSFNSLQ YFTGKIPKAT PELGSENSA SSPPRFKTEK MESKTVLPFT DKCESSNLTV NTSFGSCPQS ETIISLKID ATLTRKNHSE HPLCSEEALD PDIELSLVSE EDKQSTSNRD FETEAEDESI YFTPELYDPE DTDEEKNDLA ETDRGNRLAN NSDCILAKDL FEIRTIKEVD SAREVKAEDC IDTKLNGILH IEESKIDDID GNVKTTWINE LELGKTHEIE IKNFKPSPSK

image



Background

Alternative Names	Fanconi anemia group J protein (Protein FACJ) (EC 3.6.4.13) (ATP-dependent RNA helicase BRIP1) (BRCA1-associated C-terminal helicase 1) (BRCA1-
--------------------------	---

Product Name: Recombinant Human BRIP1 Protein
Catalog #: PEH8429



interacting protein C-terminal helicase 1) (BRCA1-interacting protein 1)

Note

For Research Use Only , Not for Diagnostic Use