


Product Name: Recombinant Human ASPSCR1 Protein
Catalog #: PEH8431



Summary

Name	ASPC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BZE9
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RPVDREPVDR EPVCHPDLE ERLQAWPAEL PDEFFELTVD DVRRRLAQLK SERKRLEEAP LVTKAFREAQ IKEKLERYPK VALRVLPDR YVLQGFFRPS ETVGDLRDFV RSHLGNPELS FYLFITPPKT VLDDHTQTLF QANLFPAALV HLGAEEPAGV YLEPGLLEHA ISPSAADVLV ARYMSRAAGS PSPLPAPDPA PKSEPAEEG ALVPPEPIPG TAQPVKRS LG KVPKWLK LPA

image



Background

Alternative Names	Tether containing UBX domain for GLUT4 (Alveolar soft part sarcoma chromosomal region candidate gene 1 protein) (Alveolar soft part sarcoma locus) (Renal papillary cell carcinoma protein 17) (UBX domain-containing protein 9)
--------------------------	--

Product Name: Recombinant Human ASPSCR1 Protein
Catalog #: PEH8431



Note

For Research Use Only , Not for Diagnostic Use