


## Summary

<b>Name</b>	SRCN1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9C0H9
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	38 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	EGVWPPPNL LSQSPKKVTA ETDfnksVDF EMPppSPPLN LHELsgPAEG ASLTPKGGNP TKGLDTPGKR SVDKAVSVEA AERDwEEKRA ALTQYSAKDI NRLLEETQAE LLKAIPDLDC ASKAHPGPAP TPDHKPPKAP HGQKAAPRTE PSGRRGSDEL TVPRYRTEKP SKSPppPPPR RSFPSSHGLT TTRTGEVVVT SKKDSAFIKK AESEELVQK PQVKLRRAVS EVARPASTPP IMASAIKDED DEDRIIAELE SGGGSVPPMK VVTPGASRLK AAQGQAGSPD KSKHGKQRAE YMRIQAQQQA TKPSKEMSGS NETSSPVSEK

## image



## Background

<b>Alternative Names</b>	SRC kinase signaling inhibitor 1 (SNAP-25-interacting protein) (SNIP) (p130Cas-associated protein) (p140Cap)
--------------------------	---

**Product Name: Recombinant Human SRCIN1 Protein**  
**Catalog #: PEH8432**

---



**Note**

For Research Use Only , Not for Diagnostic Use