


Product Name: Recombinant Human PNN Protein
Catalog #: PEH8434



Summary

Name	PININ_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H307
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QLEKAKESLK NVDENIRKLT GRDPNDVRPI QARLLALSGP GGGRGRGSLL LRRGFSDSGG GPPAKQRDLE GAVSRLGGER RTRRESRQES DPEDDDVKKP ALQSSVVATS KERTRRDLIQ DQNMDEKQKQ RNRRIFFLLM GTLQKFKQES TVATERQKRR QEIEQKLEVQ AEEERKQVEN ERRELFERR AKQTELRLL QKVELAQLQE

image



Background

Alternative Names	Pinin (140 kDa nuclear and cell adhesion-related phosphoprotein) (Desmosome-associated protein) (Domain-rich serine protein) (DRS protein) (DRSP) (Melanoma metastasis clone A protein) (Nuclear protein SDK3) (SR-like protein)
--------------------------	--

Product Name: Recombinant Human PNN Protein
Catalog #: PEH8434



Note

For Research Use Only , Not for Diagnostic Use