


Product Name: Recombinant Human ISYNA1 Protein
Catalog #: PEH8438



Summary

Name	INO1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NPH2
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QNTLVPGALE LAWQHRVFVG GDDFKSGQTK VKSVLVDFLI GSGLKTMSIV SYNHLGNNDG ENLSAPLQFR SKEVSKSNVV DDMVQSNPVL YTPGEEDHC VVIKYVPYVG DSKRALDEYT SELMLGGTNT LVLHNTCEDS LLAAPIMLDL ALLTELCQRV SFCTDMDPEP QTFHPVLSLL SFLFKAPLVP PGSPVVNALF RQRSCIENIL RACVGLPPQN

image



Background

Alternative Names	Inositol-3-phosphate synthase 1 (IPS 1) (EC 5.5.1.4) (Myo-inositol 1-phosphate synthase) (MI-1-P synthase) (MIP synthase) (hIPS) (Myo-inositol 1-phosphate synthase A1) (hINO1)
--------------------------	---

Product Name: Recombinant Human ISYNA1 Protein
Catalog #: PEH8438



Note

For Research Use Only , Not for Diagnostic Use