


Product Name: Recombinant Human BIRC6 Protein
Catalog #: PEH8439



Summary

Name	BIRC6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NR09
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NYATRIVTLE EEPIKIQHIK DPQDTITSLI LLPPDILDNR EDDCEEPIED MQLTSKNGFE REKTSDISTL GHLVITTQGG YVKILDLSNF EILAKVEPPK KEGTEEQDTF VSVIYCSGTD RLCACKGGE LHFLQIGGTC DDIDEADILV DGSLSKGIEP SSEGSKPLSN

image



Background

Alternative Names	Baculoviral IAP repeat-containing protein 6 (EC 2.3.2.27) (BIR repeat-containing ubiquitin-conjugating enzyme) (BRUCE) (RING-type E3 ubiquitin transferase BIRC6) (Ubiquitin-conjugating BIR domain enzyme apollon) (APOLLON)
--------------------------	---

Product Name: Recombinant Human BIRC6 Protein
Catalog #: PEH8439



Note

For Research Use Only , Not for Diagnostic Use