



Product Name: Recombinant Human RFOX1 Protein
Catalog #: PEH8442



Summary

Name	RFOX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NWB1
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NQEAAAAPDT MAQPYASAQF APPQNGIPAE YTAPHPHPAP EYTGQTTVPE HTLNLYPPAQ THSEQSPADT SAQTVSGTAT QTDDAAPT DG QPQTQPSENT ENKSQPKRLH VSNIPFRFRD PDLRQMFGQF GKILDVEIIF NERGSKGFGF VTFENSADAD RAREKLHGT VEGRKIEVNN ATARVMTNKK

image



Background

Alternative Names RNA binding protein fox-1 homolog 1 (Ataxin-2-binding protein 1) (Fox-1 homolog A) (Hexaribonucleotide-binding protein 1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RBFOX1 Protein
Catalog #: PEH8442

