

Product Name: Recombinant Human NDOR1 Protein
Catalog #: PEH8446



Summary

Name	NDOR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UHB4
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GSQTGTAQDV SERLGREARR RRLGCRVQAL DSYPVVNLIN EPLVIFVCAT TGQGDPPDNM KNFWRFIFRK NLPSTALCQM DFAVLGLGDS SYAKFNFAK KLHRRLLQLG GSALLPVCLG DDQHELGPDA AVDPWLRDLW

image



Background

Alternative Names NADPH-dependent diflavin oxidoreductase 1 (EC 1.18.1.-) (NADPH-dependent FMN and FAD-containing oxidoreductase) (Novel reductase 1)

Note

For Research Use Only , Not for Diagnostic Use