


**Product Name: Recombinant Human RAB23 Protein**  
**Catalog #: PEH8450**



## Summary

|                              |  |
|------------------------------|--|
| <b>Name</b>                  | RAB23_HUMAN  |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE  |
| <b>Endotoxin level</b>       | Please contact with the lab for this information   |
| <b>Tag</b>                   | C-His  |
| <b>Accession #</b>           | Q9ULC3   |
| <b>Host</b>                  | E.coli   |
| <b>Species</b>               | Human  |
| <b>Molecular Mass</b>        | 17 kDa   |
| <b>Buffer</b>                | Phosphate buffered saline  |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.   |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| <b>Sequence</b>              | IQVNDEDVRL MLWDTAGQEE FDAITKAYR GAQACVLVFS TTDRESFEAV<br>SSWREKVVAE VGDIPTVLVQ NKIDLLDDSC IKNEEAEALA KRLKLFYRT<br>SVKEDLNVNE VFKYLAEKYL QKLKQQAED  |

## image



## Background

**Alternative Names** Ras-related protein Rab-23

## Note

For Research Use Only , Not for Diagnostic Use