


Product Name: Recombinant Human MYO5B Protein
Catalog #: PEH8451



Summary

Name	MYO5B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9ULV0
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KRVRELEQER KKLQVQLEKR EQQDSKKVQA EPPQTDIDL PNADLAYNSL KRQELESENK KLNLDLNLRL KAVADQATQN NSSHGSPDSY SLLNQLKLA HEELEVRKEE VLILRTQIVS ADQRRLAGRN AEPNINARSS WPNSEKHVDQ EDAIEAYHGV CQTNSKTEDW GYLNEDGELG LAYQGLKQVA RLLEAQLQAQ SLEHEEEVEH LKAQLEALKE EMDKQQQTFC QTLLLSPEAQ VEEGVQQEIS RLTNENLDLK ELVEKLEKNE

image



Background

Alternative Names	Unconventional myosin-Vb
--------------------------	--------------------------

Product Name: Recombinant Human MYO5B Protein
Catalog #: PEH8451



Note

For Research Use Only , Not for Diagnostic Use