


Product Name: Recombinant Human MACF1 Protein
Catalog #: PEH8453



Summary

Name	MACF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UPN3
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KIEDEVTRQV AQCKCAKRFQ VEQIGENKYR FGDSQQQLRLV RILRSTVMVR VGGGWMALDE FLVKNDPCRA RGRTNIELRE KFILPEGASQ GMTPFRRSRGR RSKPSSRAAS PTRSSSSASQ SNHSCTSMPS SPATPASGTK VIPSSGSKLK RPTPTFHSSR TSLAGDTSNS SSPASTGAKT NRADPKKSAS RPSRAGSRA GSRASSRRGS DASDFDLET QSACSDTSES SAAGGQGNRS RGLNKPSKIP TMSKKTITAS

image



Background

Alternative Names	Microtubule-actin cross-linking factor 1, isoforms 1/2/3/5 (620 kDa actin-binding protein) (ABP620) (Actin cross-linking family protein 7) (Macrophin-1) (Trabeculin-alpha)
--------------------------	---

Product Name: Recombinant Human MACF1 Protein
Catalog #: PEH8453



Note

For Research Use Only , Not for Diagnostic Use