


Product Name: Recombinant Human DLEC1 Protein
Catalog #: PEH8455



Summary

Name	DLEC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y238
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EKSQDPGEL TDLTAQHFIK FEPENLRSTA RKQLIIRNAT HVELAFYWQI MKPNLQPLMP GETFSMDSIK CYPDKETAFL IMPRKGVLSP HTDHEFILSF SPHELRFHLS VLQMVLEEVP EPVSSEAESL GHSSYSVDDV IVLEIEVKGS VEPFQVLLEP YALIIPGENY IGINVKKAFK MWNNNSKSPIR YLWGKISDCH IIEVEPGTGV IEPSEVGDFE LNFTGGVPGP

image



Background

Alternative Names	Deleted in lung and esophageal cancer protein 1 (Deleted in lung cancer protein 1) (DLC-1)
--------------------------	--

Note

Product Name: Recombinant Human DLEC1 Protein
Catalog #: PEH8455



For Research Use Only , Not for Diagnostic Use