

Product Name: Recombinant Human MAP3K2 Protein
Catalog #: PEH8458



Summary

Name	M3K2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y2U5
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EFIPESMDQM LDPLSLSSPE NSGSGSCPSL DSPLDGESYP KSRMPRAQSY PDNHQEFSDY DNPIFEKFGK GGTYP RRYHV SYHHQEYNDG RKTFFPRARRT QGTSLRSPVS FSPTDHSLST

image



Background

Alternative Names Mitogen-activated protein kinase kinase kinase 2 (EC 2.7.11.25) (MAPK/ERK kinase kinase 2) (MEK kinase 2) (MEKK 2)

Note

For Research Use Only , Not for Diagnostic Use