


**Product Name: Recombinant Human ARFGEF1 Protein**  
**Catalog #: PEH8460**



## Summary

<b>Name</b>	BIG1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9Y6D6
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	44 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	SQLRKACEVA LEEIKAETEK QSPPHGEAKA GSSTLPPVKS KTNFIEADKY FLPFELACQS KCPRIVSTSL DCLQKLIAYG HLTGNAPDST TPGKKLIDRI IETICGCFQG PQTDEGVQLQ IIKALLTAVT SQHIEIHEGT VLQAVRTCYN IYLASKNLIN QTTAKATLTQ MLNVIFARME NQALQEAKQM EKERHRQHHH LLQSPVSHHE PESPQLRYLP PQTVDHISQE HEGDLDLHTN DVDKSLQDDT EPENGSDISS AENEQTEADQ ATAAETLSKN EVLYDGENHD CEEKPQDIVQ NIVEEMVNIV VGDMGEGTTI NASADGNIGT IEDGSDSENI QANGIPGTPI SVAYTPSLPD DRLSVSSNDT QESGNSSGPS

## image



## Background

**Alternative Names** Brefeldin A-inhibited guanine nucleotide-exchange protein 1 (Brefeldin A-

**Product Name: Recombinant Human ARFGEF1 Protein**  
**Catalog #: PEH8460**



---

inhibited GEP 1) (ADP-ribosylation factor guanine nucleotide-exchange factor 1) (p200 ARF guanine nucleotide exchange factor) (p200 ARF-GEP1)

**Note**

For Research Use Only , Not for Diagnostic Use