


Product Name: Recombinant Human AKAP10 Protein
Catalog #: PEH8462



Summary

Name	AKA10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43572
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ADILFCESAL FYFSEYMEKE DAVNILQFWL AADNFQSQLA AKKGQYDGQE AQNDAMILYD KYFSLQATHP LGFDDVVRLE IESNICREGG PLPNCFTTPL RQAWTTMEKV FLPGLSSNL YYKYLNDLIH SVRGDEFLGG NVSLTAPGSV GPPDESHPGS SDSSASQSSV KKASIKILKN FDEAIIVDAA SLDPESLYQR TYAGKMTFGR VSDLGQFIRE SEPEPDVRKS KGSMFSQAMK KVVQGNTDEA

image



Background

Alternative Names	A-kinase anchor protein 10, mitochondrial (AKAP-10) (Dual specificity A kinase-anchoring protein 2) (D-AKAP-2) (Protein kinase A-anchoring protein 10) (PRKA10)
--------------------------	---

Product Name: Recombinant Human AKAP10 Protein
Catalog #: PEH8462



Note

For Research Use Only , Not for Diagnostic Use