


Product Name: Recombinant Human AKAP7 Protein
Catalog #: PEH8463



Summary

Name	AKA7A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43687
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MGQLCCFPFS RDEGKISELE SSSSAVLQRY SKDIPSWSSG EKNGGEPDDA ELVRLSKRLV ENAVLKAVQQ YLEETQNKNK PGEGSSVKTE AADQNGNDNE NNRK

image



Background

Alternative Names	A-kinase anchor protein 7 isoforms alpha and beta (AKAP-7 isoforms alpha and beta) (A-kinase anchor protein 18 kDa) (AKAP 18) (Protein kinase A-anchoring protein 7 isoforms alpha/beta) (PKA7 isoforms alpha/beta)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human AKAP7 Protein
Catalog #: PEH8463

