


Product Name: Recombinant Human ACAP1 Protein
Catalog #: PEH8469



Summary

Name	ACAP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15027
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PVPPKPSIRP RPGSLRSKPE PPSEDLGSLH PGALLFRASG HPPSLPTMAD ALAHGADVNW VNGGQDNATP LIQATAANSL LACEFLLQNG ANVNQADSAG RGPLHHATIL GHTGLACLFL KRGADLGARD SEGRDPLTIA METANADIVT LLRLAKMREA EAAQGQAGDE TYLDIFRDFS LMASDDPEKL

image



Background

Alternative Names	Arf-GAP with coiled-coil, ANK repeat and PH domain-containing protein 1 (Centaurin-beta-1) (Cnt-b1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ACAP1 Protein
Catalog #: PEH8469

