


Product Name: Recombinant Human NAT2 Protein
Catalog #: PEH8472



Summary

Name	ARY2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	P11245
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LDQIRREQYI TNKEFLNSHL LPKKKHQKIY LFTLEPRTIE DFESMNTYLQ TSPTSSFITT SFCSLQTPEG VYCLVGFIIT YRKFNKYDNT

image



Background

Alternative Names	Arylamine N-acetyltransferase 2 (EC 2.3.1.5) (Arylamide acetylase 2) (N-acetyltransferase type 2) (NAT-2) (Polymorphic arylamine N-acetyltransferase) (PNAT)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use