


**Product Name: Recombinant Human NAIP Protein**  
**Catalog #: PEH8473**



## Summary

<b>Name</b>	BIRC1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q13075
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	39 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	SQVPTIDQDY ASAFEPMNEW ERNLAEKEDN VKSYMMDMQR ASPDLSTGYW KLSPKQYKIP CLEVDVNDID VVGQDMLEIL MTVFSASQRI ELHLNHSRGF IESIRPALEL SKASVTKCSI SKLELSAAEQ ELLLTLPSE SLEVSGTIQS QDQIFPNLDK FLCLKELSVD LEGNINVFSV IPEEFPNFHH MEKLLIQISA EYDPSKLVKL IQNSPNLHVF HLKCNFFSDF GSLMTMLVSC KKLTEIKFSD SFFQAVPFVA SLPNFISLKI LNLEGQQFPD EETSEKFAYI LGSLSNLEEL ILPTGDGIYR VAKLIQQCCQ

## image



## Background

<b>Alternative Names</b>	Baculoviral IAP repeat-containing protein 1 (Neuronal apoptosis inhibitory protein)
--------------------------	-------------------------------------------------------------------------------------

**Product Name: Recombinant Human NAIP Protein**  
**Catalog #: PEH8473**



---

**Note**

For Research Use Only , Not for Diagnostic Use