



Product Name: Recombinant Human BBS12 Protein
Catalog #: PEH8474



Summary

Name	BBS12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZW61
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KPQTKVEADN NTSRTLKNSL LADTCCRQSI LIHSRHFNRT DNTEGVSKPD GFQEHVTATH KTYRCNDLVE LAVGLSHGDH SSMKLVEEAV QLQYQNACVQ QGNCTKPFMF DISRIFTCCL PGLPETSSCV CPGYITVVSV SNNPVIKELQ NQPVRIVLIE GDLTENYRHL GFNKSANIKT VLDSMRLQED

image



Background

Alternative Names Bardet-Biedl syndrome 12 protein

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human BBS12 Protein
Catalog #: PEH8474

