


Summary

Name	ACSF4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q4L235
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QVAEELQQVE SCAVTWYNQE KLILFMVSKD ASVKEYIFKE LQKYLPSHAV PDELVLIDSL PFTSHGKIDV SELNKIYLNLY INLKSENKLS GKEDLWEKIQ YLWKSTLNLP EDLLRVPDES LFLNSGGDSL KSIRLLSEIE KLVGTSVPGL LEIILSSSIL EIYNHILQTV

image



Background

Alternative Names	Beta-alanine-activating enzyme (EC 6.2.1.-) (Acyl-CoA synthetase family member 4) (Protein NRPS998)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human AASDH Protein
Catalog #: PEH8475

