


Product Name: Recombinant Human BLOC1S1 Protein
Catalog #: PEH8477



Summary

Name	BL1S1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P78537
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SFRSRRGPGV PSPQPDVTML SLLKEHQAK QNERKELQEK RRREAITAAT CLTEALVDHL NVGVAQAYMN QRKLDHEVKT

image



Background

Alternative Names Biogenesis of lysosome-related organelles complex 1 subunit 1 (BLOC-1 subunit 1) (GCN5-like protein 1) (Protein RT14)

Note

For Research Use Only , Not for Diagnostic Use