


Product Name: Recombinant Human CAPN5 Protein
Catalog #: PEH8481



Summary

Name	CAN5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	O15484
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GAWTLHEDPR QNRGGGCINH KDTFFQNPQY IFEVKKPEDE VLICIQRPK RSTRREGKGE NLAIGFDIYK VEENRQYRMH SLQHKAASSI YINSRSVFLR TDQPEGRYVI IPTTFEPGHT GEFLLRVFTD VPSNCRELRL

image



Background

Alternative Names Calpain-5 (EC 3.4.22.-) (Calpain htra-3) (New calpain 3) (nCL-3)

Note

For Research Use Only , Not for Diagnostic Use