


Product Name: Recombinant Human CARD17P Protein
Catalog #: PEH8483



Summary

Name	CAR17_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5XLA6
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KQFIRSVGEG TINGLLGELL ETRVLSQEEI EIVKCENATV MDKARALLDS VIRKGAPACQ ICITYICEED SHLAGTLGLS AGPTSGNHLT

image



Background

Alternative Names Caspase recruitment domain-containing protein 17 (Caspase-1 inhibitor INCA) (Inhibitory caspase recruitment domain protein)

Note

For Research Use Only , Not for Diagnostic Use