


Product Name: Recombinant Human CNOT4 Protein
Catalog #: PEH8484



Summary

Name	CNOT4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95628
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKHQEYEQKL LQELYKLNPN FLQLSTGSVD KNKNKVTPLQ RYDTPIDKPS DSLISNGDN SQQISNSDTP SPPPGLSKSN PVIPISSNH

image



Background

Alternative Names	CCR4-NOT transcription complex subunit 4 (EC 2.3.2.27) (CCR4-associated factor 4) (E3 ubiquitin-protein ligase CNOT4) (Potential transcriptional repressor NOT4Hp) (RING-type E3 ubiquitin transferase CNOT4)
--------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Note

For Research Use Only , Not for Diagnostic Use