


Summary

Name	CD2B2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95400
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VGDEEDEDEI IVPKKKLVDP VAGSGGPGSR FKGKHS LSDS EEEDDDGGG SKYDILASED VEGQEAATLP SEGGVRITPF NLQEEMEEGH FDADGNYFLN RDAQIRDSWL DNIDWVKIRE RPPGQRQASD SEEDSLGQT SMSAQALLEG LLELLPRET VAGALRRLGA RGGGKGRKGP GQPSSPQRLD RLSGLADQMV ARGNLGVYQE TRERLAMRLK GLGCQTLGPH NPTPPPSLDM FAEELAEEL ETPTPTQRGE

image



Background

Alternative Names	CD2 antigen cytoplasmic tail-binding protein 2 (CD2 cytoplasmic domain-binding protein 2) (CD2 tail-binding protein 2) (U5 snRNP 52K protein) (U5-52K)
--------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------

Product Name: Recombinant Human CD2BP2 Protein
Catalog #: PEH8485



Note

For Research Use Only , Not for Diagnostic Use