


Product Name: Recombinant Human CENPU Protein
Catalog #: PEH8488



Summary

Name	CENPU_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q71F23
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SGNEASEIES VKISAKKPGR KLRPISDDSE SIEESDTRRK VKSAEKISTQ RHEVIRTAS SELSEKPAES VTSKKTGPLS AQPSVEKENL AIESQSKTQK KKGISHDKRK KSRSKAIGSD

image



Background

Alternative Names	Centromere protein U (CENP-U) (Centromere protein of 50 kDa) (CENP-50) (Interphase centromere complex protein 24) (KSHV latent nuclear antigen-interacting protein 1) (MLF1-interacting protein) (Polo-box-interacting protein 1)
--------------------------	---

Note

Product Name: Recombinant Human CENPU Protein
Catalog #: PEH8488



For Research Use Only , Not for Diagnostic Use