


**Product Name: Recombinant Human CHD3 Protein**  
**Catalog #: PEH8490**



## Summary

<b>Name</b>	CHD3_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q12873
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	17 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	ALPPPPAADI QPPPIRRAKT KEGKGPCHKR RSKSPRVPDG RKKLRGKKMA PLKIKLGLLG GKRKKGGSYV FQSDEGPEPE AEESDLDSGS VHSASGRPDG PVRTKCLKRG RPGRKCKKVL GCPAVAGEEE VDGYETDHQD

## image



## Background

**Alternative Names** Chromodomain-helicase-DNA-binding protein 3 (CHD-3) (EC 3.6.4.12) (ATP-dependent helicase CHD3) (Mi-2 autoantigen 240 kDa protein) (Mi2-alpha) (Zinc finger helicase) (hZFH)

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human CHD3 Protein**  
**Catalog #: PEH8490**

---

