


**Product Name: Recombinant Human KIF4A Protein**  
**Catalog #: PEH8491**



## Summary

<b>Name</b>	KIF4A_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O95239
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	43 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	VETLEDQELK ENVEIICNLQ QLITQLSDET VACMAAAIDT AVEQEAQVET SPETSRSSDA FTTQHALRQA QMSKELVELN KALALKEALA RKMTQNDSQL QPIQYQYQDN IKELELEVIN LQKEKEELVL ELQTAKKDAN QAKLSERRRK RLQELEGQIA DLKKKLNEQS KLLKLKESTE RTVSKLNQEI RMMKNQRVQL MRQMKEDA EK FRQWKQKKDK EVIQLKERDR KRQYELLKLE RNFQKQSNVL RRKTEEAAAA NKRLKDALQK QREVADKRKE TQSRGMEGTA ARVKNWLGNE IEVMVSTEEA KRHLNDLLED RKILAQDVAQ LKEKKESGEN PPPKLRRTTF

## image



## Background

**Alternative Names** Chromosome-associated kinesin KIF4A (Chromokinesin-A)

**Product Name: Recombinant Human KIF4A Protein**  
**Catalog #: PEH8491**

---



**Note**

For Research Use Only , Not for Diagnostic Use