


Product Name: Recombinant Human COG2 Protein
Catalog #: PEH8494



Summary

Name	COG2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14746
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EFMKEDFDVD HFVSDCRKRV QLEELRDDLE LYYKLLKTAM VELINKDYAD FVNLSTNLVG MDKALNQLSV PLGQLREEVL SLRSSVSEGI RAVDERMSKQ

image



Background

Alternative Names	Conserved oligomeric Golgi complex subunit 2 (COG complex subunit 2) (Component of oligomeric Golgi complex 2) (Low density lipoprotein receptor defect C-complementing protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use