


Product Name: Recombinant Human CUL7 Protein
Catalog #: PEH8495



Summary

Name	CUL7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14999
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VEADEYQGAV ASRVLGRALP AWRWRPMTTEL YAVPYVLPED EDTEECEHLT LAEWWELEFF IKKLDGPDHQ EVLQILQENL DGEILDDEIL AELAVPIELA QDLLLLTLPQR LNSALRDLI NCHVYKKGYP EALAGNQAYP SLLEAQEDVL LLDAQAQAKD SEDAQKVEAK EPPSQSPNTP LQRLVEGYGP AGKILLDLEQ ALSSEGTQEN KVKPLLLQLQ RQPQPFLALM QSLDTPETNR TLHLTVLRIL

image



Background

Alternative Names Cullin-7 (CUL-7)

Note

Product Name: Recombinant Human CUL7 Protein
Catalog #: PEH8495



For Research Use Only , Not for Diagnostic Use