


Product Name: Recombinant Human HCCS Protein
Catalog #: PEH8498



Summary

Name	CCHL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P53701
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PSGPTCEKKT YSVPAHQERA YEYVECPiRG TAAENKENLD PSNLMPPPNQ TPAPDQPFAL STVREESSIP RADSEKKWVY PSEQMFWNAM LKKGWKWKDE DISQKDMYNI IRIHNQNNEQ AWKEILKWEA LHAAECPCGP SLIRFGGKAK EYSPRARIRS WMGYELPFDR HDWIINRCGT EVRYVIDYYD GGEVNKDYQF TILDVRPALD

image



Background

Alternative Names	Cytochrome c-type heme lyase (CCHL) (EC 4.4.1.17) (Holocytochrome c-type synthase)
--------------------------	--

Note

Product Name: Recombinant Human HCCS Protein
Catalog #: PEH8498



For Research Use Only , Not for Diagnostic Use