


Summary

Name	CP4AB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q02928
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KAQLQKEGEL EKIKRKRHLD FLDILLAKM ENGSILSDKD LRAEVDTFMF EGHDTTASGI SWILYALATH PKHQERCREE IHSLLGDGAS ITWNHLDQMP YTTMCIKEAL RLYPPVPGIG RELSTPVTFP

image



Background

Alternative Names	Cytochrome P450 4A11 (EC 1.14.14.1) (20-hydroxyeicosatetraenoic acid synthase) (20-HETE synthase) (CYP4A11) (CYP1VA11) (Cytochrome P-450HK-omega) (Cytochrome P450HL-omega) (Fatty acid omega-hydroxylase) (Lauric acid omega-hydroxylase) (Long-chain fatty acid omega-monooxygenase) (EC 1.14.14.80)
--------------------------	--

Product Name: Recombinant Human CYP4A11 Protein
Catalog #: PEH8499



Note

For Research Use Only , Not for Diagnostic Use