


Product Name: Recombinant Human CYTH3 Protein
Catalog #: PEH8500



Summary

Name	CYH3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43739
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GVPEDLSLEE REELDIRRR KKEIIDIER LKYEIAEVMT EIDNLTVEE SKTTQRNKQI AMGRKKFNMD PKKGIQFLIE NDLLQSSPED VAQFLYKGEG LNKTVIGDYL GERDEFNIKV LQAFVELHEF ADLNLVQALR QFLWSFRLPG EAQKIDRMME AFASRYCLCN PGVFQSTDTC YVLSFAIIML NTSLNHNHNR DKPTAERFIA MNRGINEGGD LPEELLRNLY ESIKNEPFKI PEDDGNDLTH

image



Background

Alternative Names	Cytohesin-3 (ARF nucleotide-binding site opener 3) (Protein ARNO3) (General receptor of phosphoinositides 1) (Grp1) (PH, SEC7 and coiled-coil domain-containing protein 3)
--------------------------	--

Product Name: Recombinant Human CYTH3 Protein
Catalog #: PEH8500



Note

For Research Use Only , Not for Diagnostic Use