


Product Name: Recombinant Human CLNK Protein
Catalog #: PEH8501



Summary

Name	CLNK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z7G1
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GQPTWNTQTR LERVDKPISK DVRSQNIKGD ASVRKNKIPL PPPRPLITLP KKYQPLPPEP ESSRPPLSQR HTFPEVQRMP SQISLRDLSE VLEAEKVPHN QRKPESTHLL ENQNTQEIPL AISSSSFTTS NHSVQNRDHR GGMQPCSPQR CQPPASCSPH

image



Background

Alternative Names	Cytokine-dependent hematopoietic cell linker (Mast cell immunoreceptor signal transducer)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CLNK Protein
Catalog #: PEH8501

