


Product Name: Recombinant Human POLR3G Protein
Catalog #: PEH8505



Summary

Name	RPC7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15318
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GEGEYMLAL KQELRETMKR MPYFIETPEE RQDIERYSKR YMKVYKEEWI PDWRRLPREM MPRNKCKKAG PKPKKAKDAG KGTPLTNTED VLKKMEELEK RGDGEKSDEE NEEKEGSKEK SKEGDDDDDD DAAEQEEYDE

image



Background

Alternative Names	DNA-directed RNA polymerase III subunit RPC7 (RNA polymerase III subunit C7) (DNA-directed RNA polymerase III subunit G) (RNA polymerase III 32 kDa apha subunit) (RPC32-alpha) (RNA polymerase III 32 kDa subunit) (RPC32)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human POLR3G Protein
Catalog #: PEH8505

