


**Product Name: Recombinant Human MGRN1 Protein**  
**Catalog #: PEH8506**



## Summary

|                              |                                                                                                                                                                                                                                                                                                               |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | MGRN1_HUMAN                                                                                                                                                                                                                                                                                                   |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                                           |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                                                                                                              |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                                                                                                                         |
| <b>Accession #</b>           | O60291                                                                                                                                                                                                                                                                                                        |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                                                                                                                        |
| <b>Species</b>               | Human                                                                                                                                                                                                                                                                                                         |
| <b>Molecular Mass</b>        | 31 kDa                                                                                                                                                                                                                                                                                                        |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                                                                                                                     |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                                                                            |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                                                                                |
| <b>Sequence</b>              | GVEDIDIQAN SAYRYPPKSG NYFASHFFMG GEKFDTPHPE GYLFGENMDL<br>NFLGSRPVQF PYVTPAPHEP VKTLRSLVNI RKDSLRLVRY KDDADSPTED<br>GDKPRVLYSL EFTFDADARV AITIYCQASE EFLNGRAVYS PKSPSLQSET<br>VHYKRGVSQQ FSLPSFKIDF SEWKDDELNF DLDRGVFPVV IQAVVDEGDV<br>VEVTGHAHVL LAAFEKHMDG SFSVKPLKQK QIVDRVSYLL QEIYG IENKN<br>NQETKPSDDE |

## image



## Background

|                          |                                                                                                                                                        |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | E3 ubiquitin-protein ligase MGRN1 (EC 2.3.2.27) (Mahogunin RING finger protein 1) (RING finger protein 156) (RING-type E3 ubiquitin transferase MGRN1) |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|

**Product Name: Recombinant Human MGRN1 Protein**  
**Catalog #: PEH8506**

---



**Note**

For Research Use Only , Not for Diagnostic Use