


**Product Name: Recombinant Human EVI5 Protein**  
**Catalog #: PEH8509**



## Summary

<b>Name</b>	EV15_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O60447
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	39 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	IEKRNSLPD ENNIARLQEE LIAVKLREAE AIMGLKELRQ QVKDLEEHWQ RHLARTTGRW KDPPKKNAMN ELQDELMTIR LREAETQAEI REIKQRMMEM ETQNQINSNH LRRAEQEVIS LQEKVQYLSA QNKGLLTQLS EAKRKQAEIE CKNKKEEVM AV RLREADSIAA VAELRQHIAE LEIQKEEGKL QGQLNKSDSN QYIGELKDQI AELNHELRLCL KGQRGFSGQP PFDGIHIVNH LIGDDES FHS SDEDFIDNSL QETGVGFPLH GKSGSMSLDP AVADGSESET EDSVLETRES NQVVQKERPP RRRESYSTTV

## image



## Background

**Alternative Names** Ecotropic viral integration site 5 protein homolog (EVI-5) (Neuroblastoma stage 4S gene protein)

**Product Name: Recombinant Human EVI5 Protein**  
**Catalog #: PEH8509**

---



**Note**

For Research Use Only , Not for Diagnostic Use