

Product Name: Recombinant Human EFHC1 Protein
Catalog #: PEH8510



Summary

Name	EFHC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5JVL4
Host	E.coli
Species	Human
Molecular Mass	62 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FLPGTSFKDS TKTAFHRSQT LSYRNGYAIV RRPTVGIGGD RLQFNQLSQA ELDELASKAP VLTYGQPKQA PPADFIPAHV AFDKKVLKFD AYFQEDVPMS TEEQYRIRQV NIYYLEDDS MSVIEPVVEN SGILQGKLIK RQRLAKNDRG DHYHWKDLNR GINITYGKT FRVVDQDQFT QVFLESQGIE LNPPEKMALD PYTELKQPL RKYVTPSDFD QLKQFLTFDK QVLRFYAIWD DTDSMYGECR TYIIHYLMD DTVEIREVHE RNDGRDPFPL LMNRQVRPKV LVENAKNFPQ CVLEISDQEV LEWYTAKDFI VGKSLTILGR TFFIYDCDPF TRRYYKEKFG ITDLPRIDVS KREPPPVKQE LPPYNGFGLV EDSAQNCFAL IPKAPKKDVI KMLVNDNKVL RYLAVLESPI PEDKDRRFVF SYFLATDMIS IFEPVVRNSG IIGGKYLGR KVKPYSTVD NPVYGPSDF FIGAVIEVFG HRFIILDTDE YVLKYMESNA AQYSPEALAS

image



Background

Product Name: Recombinant Human EFHC1 Protein
Catalog #: PEH8510



Alternative Names

EF-hand domain-containing protein 1 (Myoclonin-1)

Note

For Research Use Only , Not for Diagnostic Use