


Product Name: Recombinant Human EFHC1 Protein
Catalog #: PEH8511



Summary

Name	EFHC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5JVL4
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PYTELRKQPL RKYVTPSDFD QLKQFLTFDK QVLRFYAIWD DTDSMYGECR TYIIHYLMD DTVEIREVHE RNDGRDPFPL LMNRQVRPKV LVENAKNFPQ CVLEISDQEV LEWYTAKDFI VGKSLTILGR TFFIYDCDPF TRRYKKEKFG ITDLPRIDVS KREPPPVKQE LPPYNGFGLV EDSAQNCFAL IPKAPKKDVI KMLVNDNKVL RYLAVLESPI PEDKDRRFVF SYFLATDMIS IFEPPVRNSG IIGGKYLGR T KVVKPYSTVD NPVYGPSDF FIGAVIEVFG HRFIILDTE YVLKYMESNA AQYSPEALAS IQNHVRKREA

image



Background

Alternative Names EF-hand domain-containing protein 1 (Myoclonin-1)

Product Name: Recombinant Human EFHC1 Protein
Catalog #: PEH8511



Note

For Research Use Only , Not for Diagnostic Use