


Product Name: Recombinant Human TSFM Protein
Catalog #: PEH8515



Summary

Name	EFTS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P43897
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MKLRRKTGYS FVNCKKALET CGGDLKQAEI WLHKEAQKEG WSKAAKLQGR KTKEGLIGLL QEGNTTVLVE VNCETDFVSR NLKFQLLVQQ VALGTMMHCQ TLKDQPSAYS KGFLNSEELS GLPAGPDREG SLKDQLALAI GKLGENMILK RAAWVKVPSG FYVGSYVHGA MQSPSLHKLV LGKYGALVIC ETSEQKTNLE DVGRRLGQHV VGMAPLSVGS LDDEPGGEAE TKMLSQPYLL DPSITLGQYV QPQGVSVDVDF VRFECGEGEE

image



Background

Alternative Names Elongation factor Ts, mitochondrial (EF-Ts) (EF-TsMt)

Product Name: Recombinant Human TSFM Protein
Catalog #: PEH8515



Note

For Research Use Only , Not for Diagnostic Use