



Product Name: Recombinant Human MAP7 Protein
Catalog #: PEH8517



Summary

Name	MAP7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14244
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PAMVSAPSST VNASASVKTS AGTTDPEEAT RLLAEKRRLA REQREKEERE RREQEELERQ KRELAQRVA EERTTRREEE SRRLEAEQAR EKEEQLQRQA EERALRERE AERAQRQKEE EARVREEAER VRQEREKHFQ REEQERLERK KRLEEIMKRT RRTEATDKKT SDQRNGDIAK GALTGGTEVS ALPCTTNAPG NGKPVGSPHV VTSHQSKVTV ESTPDLEKQP NENGVSVQNE NFEEIINLPI GSKPSRLDVT NSESPEIPLN PILAFDDEGT

image



Background

Alternative Names Ensconsin (Epithelial microtubule-associated protein of 115 kDa) (E-MAP-115) (Microtubule-associated protein 7) (MAP-7)

Product Name: Recombinant Human MAP7 Protein
Catalog #: PEH8517



Note

For Research Use Only , Not for Diagnostic Use